**ASPIRATORE VEGA UNO  
ASPIRATOR VEGA UNO  
ASPIRATEUR VEGA UNO  
ABSAUGER VEGA UNO  
ASPIRADOR VEGA UNO**

**MANUALE D'USO  
USER MANUAL  
MODE D'EMPLOI  
HANDBUCH  
MANUAL DE INSTRUCCIONES**

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**VEGA UNO** is an electric device for medical use working with 230V ~ / 50 Hz network electricity, to be used for the nasal aspiration, oral aspiration, tracheal aspiration of the body liquids (mucus or catarrh) in the adult or in the children. Easily transportable from one hospital ward to another, particularly suitable for tracheotomized patients, minor surgical applications and post-operative therapy at home.  
Easily portable equipment designed for not continuous use 20 min. ON / 40 min. OFF.  
Made of highly heat-resistant, electrically insulated plastic material in conformity with the latest European safety standards.  
Supplied with a complete polycarbonate autoclavable jar complete with overflow valve.  
Equipped with aspiration regulator and vacuum indicator located on the front panel.

### GENERAL WARNING

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE**

**ONLY HIGHLY QUALIFIED STAFF USE RESERVED**

**THE INSTRUMENT MUST NOT BE DISASSEMBLED  
FOR A TECHNICAL SERVICE ALWAYS CONTACT GIMA S.p.A.**


### IMPORTANT SAFETY RULES

1. On opening the packaging, check the integrity of the appliance, paying particular attention to the presence of damage to the plastic parts, which may make access possible to internal live parts and also to breakage and / or peeling of the power supply cable. **In these cases don't connect the plug to the electric socket.**  
**Carry out these controls before each use;**
2. before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected;
3. If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and / or extension adapters is not generally recommended. Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions;
4. Respect the safety regulations indicated for electrical appliances and particularly:
  - Use original components and accessories provided by the manufacturer to guarantee the highest efficiency and safety of the device;
  - The device can be used only with the bacteriological filter;
  - Never immerse the appliance into water;
  - Position the appliance on flat stable surfaces;
  - Position the device in a way that the air inlets on the back aren't obstructed;
  - Never use the device in environments which have anaesthetic mixtures inflammable with air, oxygen or nitric oxide;
  - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
  - Keep off the reach of children or not capable people without supervision;
  - Don't leave the appliance connected to the power supply socket when not in use;
  - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
  - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources;
  - Don't use the device thoracic drainage.
5. For repairs, exclusively contact GIMA technical service and request the use of original spare parts.  
Failure to comply with the above can jeopardise the safety of the device;
6. **This medical device must be destined exclusively for the use for which it has been designed and described in this manual.** Any different use must be considered incorrect and therefore dangerous; the manufacturer cannot be considered liable for damage caused by improper, incorrect and / or unreasonable use or if the appliance is used in electrical plants that are not in compliance with the regulations in force;
7. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents;
8. Instrument and accessories discharging must be done following current law regulations in every country of use.
9. None of electric or mechanical parts have been designed to be repaired by customers or end-users. Don't open the device, do not mishandle the electric / mechanical parts. Always contact GIMA technical assistance
10. Using the device in environmental conditions different than those indicated in this manual may harm seriously the safety and the technical characteristics of the same.

## TECHNICAL CHARACTERISTICS

TYPOLOGY (MDD 93/42/EEC)	Dispositivo Medico Classe IIa	
MODEL	<b>VEGA UNO</b>	
UNI EN ISO 10079-1	HIGH VACUUM / LOW FLOW	
POWER FEEDING	230V~ / 50Hz	110V~ / 60Hz (no CE 0123)
POWER CONSUMPTION	184VA	105VA
FUSE	F 1 x 1.6A 250V	F 1 x 4A 250V
MAXIMUM SUCTION PRESSURE (without jar)	-75kPa (-0.75 Bar) Regolable from -75kPa (-0.75bar) to -10kPa (-0.10 bar)	
MAXIMUM SUCTION FLOW (without jar)	15 l/min	
WEIGHT	2.2 Kg	
SIZE	255 x 190 (h) x 165 mm	
DUTY CYCLE (to 35°C and 110% operating voltage)	20 min ON / 40 min OFF	
SICILICONE TUBE SIZE	Ø 6 x 10 mm	
ACCURACY OF VACUUM INDICATOR	± 5%	
WORKING CONDITION	Room temperature:	10 ÷ 40°C
	Room humidity percentage:	20 ÷ 85% RH
	Altitude:	0 ÷ 2000m s.l.m.
CONSERVATION CONDITION AND TRASPORT	Room temperature:	-40÷ 70°C
	Room humidity percentage:	10 ÷ 95% RH

## SYMBOLS

	Class II isolation equipment
<b>CE 0123</b>	CE marking in conformity with EC directive 93/42/EEC Manufactured by: CA.MI. snc di Attolini Mario & C. Via Ugo La Malfa nr.31 – 43010 Pilastrò (PR) Italia
	Warning, consult the instruction manual
	To Preserve in place coolness and dry land
	Conservation temperature: -40 ÷ 70°C
	Type B equipment
	Fuse
~	Alternate Current
Hz	Mains Frequency
I	ON
O	OFF

<b>Guidance and manufacturer's declaration – Electromagnetic Emissions</b>		
The surgical aspirator <b>VEGA UNO</b> is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator <b>VEGA UNO</b> should assure that it's used in such an environment.		
<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
Power disturbance CISPR11	Group 1	The surgical aspirator VEGA UNO only used RF energy only for its internal functioning. Therefore its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The surgical aspirator VEGA UNO can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes.
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

<b>Guidance and manufacturer's declaration – Electromagnetic Immunity</b>		
The surgical aspirator <b>VEGA UNO</b> is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator <b>VEGA UNO</b> should assure that it's used in such an environment.		
<b>Immunity Test</b>	<b>Compliance</b>	<b>Electromagnetic environments - guidance</b>
Electrostatic discharge (ESD) IEC/EN 61000-4-2	± 6kV on contact ± 8kV in air	Floors should be wood, concret or ceramic tile. If floors are coverei with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	± 2kV power supply	Mains power quality should be that of a typical commercial environment or hospital
Surge IEC/EN 61000-4-5	± 1kV differential mode	Mains power quality should be that of a typical commercial environment or hospital
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	5%U <sub>T</sub> for 0.5 cycle 40%U <sub>T</sub> for 05 cycle 70%U <sub>T</sub> for 25 cycle <5%U <sub>T</sub> for 5 sec	Mains power quality should be that of a typical commercial environment or hospital If the user of the surgical aspirator VEGA UNO. request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Conducted Immunity IEC/EN 61000-4-6	3Vrms 150kHz to 80MHz (for appliances that aren't life - supporting)	-
Irradiated Conducted IEC/EN 61000-4-3	3V/m 80MHz to 2.5 GHz (for appliances that aren't life - equipment)	-
Note U <sub>T</sub> is the value of the power supply voltage		

## ACCESSORIES SUPPLIED

DESCRIPTION
COMPLETE ASPIRATION JAR 1000cc
CONICAL FITTING
TUBES SET 6 mm x 10 mm
ASPIRATION PROBE CH20
ANTIBACTERIAL FILTER

The filter is produced with (PTFE) hydrophobic material which prevents fluids entering the pneumatic circuit. When the filter is wet, it's not possible to use the unit therefore the filter should be changed immediately.

In case of possible contamination or discolouration, change the filter immediately.

Don't use the suction unit without the protection filter fitted. If the suction unit is used in an emergency or in a patient where the risk of contamination is not known the filter must be changed after each use.

Available under request with different versions with complete jar 2000cc.

The suction probe must be changed at the end of each application.

## CLEANING OF ACCESSORIES

To clean the plastic housing of the device wear disposable latex gloves and clean with denaturated alcohol or hypochlorite solutions.

Washing and / or cleaning the autoclavable jar as to be carried out as follows:

- Wear protection gloves and apron (glasses and face mask if necessary) to avoid contact with contaminating substances;
- Disconnect the jar from the device
- Disconnect all tubes from the jar and the protection filter
- Empty and dispose of the content and of the suction catheter according to the laws in force in your country;
- Separate all parts of the cover (overflow valve, o-ring);

After disposing of disposable parts and disassembling the jar wash in running cold water and rinse thoroughly.

Then soak in warm water (temperature shall not exceed 60°C).

Wash thoroughly and if necessary use a non-abrasive brush to remove incrustations.

Rinse in running warm water and dry all parts with a soft cloth (non-abrasive).

The jar and the cover can be autoclaved by placing the parts into the autoclave and running one sterilization stem cycle at 121°C (1 bar relative pressure) making sure that the jar is positioned upside down.

Mechanical resistance of the jar is guaranteed up to 30 cycles of sterilization and cleaning at the indicated conditions (EN ISO 10079-1). Beyond this limit the physical-mechanical characteristics of the plastic may decrease and replacement of the part is therefore recommended.

After sterilization and cooling at environment temperature of the parts make sure that these are not damaged.

Assemble the jar as follows:

- Place the overflow valve into its seat in the cover (under VACUUM connector)
- Insert floating valve keeping the o-ring towards the opening of the cage
- Place the o-ring into its seat around the cover
- After completing assembling operations always make sure that cover seals perfectly to avoid vacuum leakages or liquid exit

The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 120°C.

The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C.

The device is ready for a new employment now.

**DO NOT WASH, STERILIZE OR PUT IN AUTOCLAVE THE ANTIBACTERIAL FILTER**

## MAINTENANCE

The **VEGA UNO** suction equipment does not need maintenance or lubrication.

It is necessary to check functioning and instrument before every use.

Unpack the instrument and **always check** integrity of plastic parts and feeding cable, they might have been damaged during previous use.

Connect cable to electrical network and turn switch on.

Close the aspiration outlet with your finger and with suction regulator in maximum vacuum position check that the vacuum indicators reaches -75 kPa (-0.75 bar) maximum.

Rotate the knob from right to left and check the aspiration regulating control.

The vacuum indicator should go down -25 kPa (-0.25 bar).

Verify that loud noises are not present, these can indicate wrong functioning.

A protection fuse (**F 1x1.6A 250V and F 1 x 4A 250V for voltage 110V/60Hz**) reachable from exterior and it situated in the plug protects the instrument. For use replacing, always check the type and the range indicated.

Fault type	Cause	Solution
1.The suction unit doesn't work	Cable is damaged  External power source failure	Replace the cable  Check the external power source
2.No aspiration	Jar Cap badly screwed down	Unscrewed the cap, then rescrew it correctly
3. No aspiration	Lid seal not in its seat	Unscrew the cap and insert the seal properly in its seat
4.The Vacuum power on the patient side is either very low or absent	a) Vacuum regulator set to minimum b) Protection filter blocked or damaged c) Connection tubes blocked, kinked or disconnected d) Shut-off valve blocked or damaged e) Pump motor damaged	a) Turn the vacuum regulator clockwise and check the value of the vacuum on the gauge b) Replace the filter c) Replace or reconnect the tubes, check the jar connections d) Empty the jar, or disconnect the tube from the jar and unblock the shut-off valve. The unit will only work in the upright position e) Refer to authorised service personnel
5. The float doesn't close	If the cap has been washed, ensure that the float is not partially detached	Insert the float into it's place
6. The float doesn't close	The float it's covered by dirty material	Unscrewed the cap, leave the and put in on autoclave
7. Low suction	Foam inside the jar	Fill the jar to 1/3 full of ordinary water
<b>Faults 1 - 2 - 3 - 4 - 5 - 6 - 7</b>	<b>None of the remedies has achieved the desired results</b>	<b>Contact the seller or GIMA S.p.A. After-sales Assistance Service</b>

If the overflow security system it's activated, don't proceed with the liquid aspiration.

If the overflow security system doesn't work there are two cases:

1° case – If the overflow security system doesn't work the aspiration will be stopped by the bacteriological filter who avoid the liquid penetration inside the device.

2° case – If both the security system doesn't work, there is the possibility that liquid comes inside the device, in this case return the device to GIMA technical service.

**The manufacturer will provide upon request electric diagrams, components list, descriptions, setting instructions and any other information that can help the technical assistance staff for product repair.**

**BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT GIMA TECHNICAL SERVICE.**

**GIMA S.p.A. DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.**

## INSTRUCTIONS

- Connect the short silicon tube, with antibacterial filter, to the suction connector.
- The other tube already connected to the filter has to be connected to the "VACUUM" jar outlet, where has been fixed the red float (security float).

When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to avoid liquid penetration inside the device.

**The device must be used on a plan of horizontal operation.**

**WARNING:** Ensure that the FLUID SIDE or IN marker on the filter is on the side facing the collection jar lid and fitted into the "VACUUM". A wrong connection causes immediate destruction in case of contact with sucked liquids.

- Connect the longer silicone tube to the "PATIENT" jar outlet
- Connect the other end of the long silicon tube to the probe plastic connector then connect the suction probe to it.
- Connect the power cord to the device then connect the plug to the electrical mains supply.
- Push switch on position I to start suction.
- Unscrew the jar's lid and fill the jar 1/3 full of ordinary water (this for an easy cleaning operations and an rapid reaching of the functionality vacuum) then rescrew the lid on the jar correctly.
- During operation the jar has to be in vertical position to avoid overflow valve to cut off aspiration. Should this happen, switch off the device and disconnect the tube from the jar cover (from "VACUUM" outlet).
- Once finished push switch on **O** position and unplug.

**NEVER USE THE DEVICE WITHOUT JAR AND / OR PROTECTION FILTER**

**MAKE SURE THAT CHILDREN AND/OR MENTALLY ILL PEOPLE DO NOT USE THE DEVICE WITHOUT ADULT SURVEILLANCE**

**ALWAYS PLACE THE DEVICE IN POSITIONS FOR EASY DISCONNECTION.**

## RULES FOR RETURNING AND REPAIRING

**COMPLYING WITH THE NEW EUROPEAN RULES, GIMA S.p.A. INDICATES THE IMPORTANT POINTS TO PROTECT INSTRUMENT AND OPERATORS HYGIENE.**

**THESE RULES MUST BE RESPECTED IN ORDER TO GUARANTEE HYGIENE AND SAFETY TO ALL THE PEOPLE OPERATING WITH THE INSTRUMENT TO OBTAIN QUALITY AND WELL BEING.**

GIMA S.p.A. warrants its products for **24 months** after purchasing date.

In front of this warranty, GIMA S.p.A. will be obliged only to repair or substitute free of charge the products or parts of them that, after verification effected on our factory, or our authorized Service Center, by the Technical Service, results defective. The product must be accompanied by a description of the defect.

The warranty, with exclusion of responsibility for direct and indirect damages, it is thought limited to the solos defects of material or workmanship and it stops having effect when the device results however gotten off, tampered or sheltered out of the Factory or from the Authorized Service center.

The commodity always travels to risk and danger of the buyer, without any responsibility of GIMA for damages caused by the transport or dismay from the vector.

Every returned instrument will be hygienically checked before repairing.

If GIMA finds instrument not suitable for repairing due to clear signs of internal or external contamination, the same will be returned to customer with specification of NOT REPAIRED INSTRUMENT, accompanied by an explanation letter.

GIMA will decide if contamination is due to bad functioning or misuse. If contamination is due to bad functioning, GIMA will substitute the instrument, only if a SALE RECEIPT and STAMPED GUARANTEE accompany the same.

GIMA is not responsible for contaminated accessories, they will be substitute at customer's expenses.

For this reason it is **COMPULSORY** to carefully disinfect the external part of the instrument and accessories with a cloth soaked in methylated spirits or hypochlorite-based solutions. Put the instrument and accessories in a bag with indication of disinfecting. We also request to specify the kind of fault, in order to speed up repairing procedures.

To this end, please read the instructions carefully in order to avoid damaging the equipment through improper use.

Always specify the fault encountered so that GIMA S.p.A. can establish whether it falls into the category of the faults covered by the guarantee.

**The manufacturer cannot be held liable for accidental or indirect damages should the device be modified, repaired without authorization or should any of its component be damaged due to accident or misuse.**

**Any minimal modification / repair on the device voids the warranty and does not guarantee the compliance with the technical requirements provided by the MDD 93/42/EEC Directive and its normatives.**

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## IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:

In respect of art. 13 Decreto Legislativo 25 Luglio 2005, n.151 "Actuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal" The symbol as over applied on the device or its packaging means that at the end of its useful life the product must not be disposed of with domestic waste. At the end of device useful, the user will must deliver it to the able collecting centres for electric and electronic garbage, or give back to the retailer in the moment of equivalent new device purchasing, one against one. Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal. It also allows the recovery of materials of witch it's made up in order to obtain an important saving of energy and resources and to avoid negative effects to the ambient and health. In case of abusive disposal of device by user, will be applied administrative endorsements in compliance with current standard.

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