
Features

- Seamless stainless steel sharp cutting edge
- Ribbed handle for sure grip and control
- Size embossed on each punch for easy identification
- Wide range of sizes
- Sterile and ready to use
- Disposable

Clinical areas

- Dermatology
- Gynecology
- Podiatry







Long type	
BP-L15K	1.5mm
BP-L20K	2mm
BP-L30K	3mm
BP-L35K	3.5mm
BP-L40K	4mm
BP-L50K	5mm
BP-L60K	6mm
BP-L80K	8mm

Package

- 20 punches in a dispenser box
- 50 boxes per shipping carton

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Features

- Stainless steel semi sharp cutting edge for soft or weak tissue
- Size embossed on each curette for easy identification
- Wide range of sizes
- Sterile and ready to use
- Disposable

Used in the following minor skin surgical procedures

- Curettage of a basal cell carcinoma
- Curettage of a verruca
- Curettage of other benign neoplasms of the skin

Clinical areas

- Dermatology
- Surgery
- Podiatry
- Veterinary

Package

- 20 curettes in a dispenser box
- 20 boxes per shipping carton

REF	SIZE	
MK402	2 mm	 Actual size
MK403	3 mm	
MK404	4 mm	
MK405	5 mm	
MK407	7 mm	

A dermal curette biopsy technique for various dermatological diseases

A punch (Trepan) biopsy technique for various dermatological diseases

References : Manabu MAEDA, MD & Ph D *Department of Dermatology, Prefectural Gifu Hospital, 4-6-1 Noishiki, Gifu City, 500-8717, Japan*

Product specifications may be changed without notice.



kai industries co., ltd.
1110, Oyana, Seki City, Gifu Pref. 501-3992 Japan
Phone: +81 (0) 575 28 6600 Fax: +81 (0) 575 28 6611
Email: info-ind@kai-group.com http://www.kai-ind.co.jp

kai Europe GmbH
Löhndorfer Str. 51, D-42699 Solingen, Germany
Phone: +49 (0) 212 232380 Fax: +49 (0) 212 2323899

DISPOSABLE BIOPSY PUNCH WITH PLUNGER SYSTEM

STERILE • STAINLESS STEEL • SINGLE USE • COLOR CODED • PATENT PENDING

	BPP-10F 1.0mm w/plunger	
	BPP-15F 1.5mm w/plunger	
	BPP-20F 2.0mm w/plunger	
NEW	BPP-30F 3.0mm w/plunger	
NEW	BPP-40F 4.0mm w/plunger	

Punches are available in three sizes, marked clearly on the handle, blister pack and box
Punches are individually packed in color-coded, transparent, sturdy blister packs, 20 per dispenser box

Product specifications may be changed without notice

kai industries co.,ltd.
1110, Oyana, Seki City, Gifu Pref. 501-3992 Japan
Phone : +81 (0) 575 28 6600 Fax : +81 (0) 575 28 6611
Email: info-ind@kai-group.com http://www.kai-ind.co.jp

kai Europe GmbH
Löhdorfer Str. 51, D-42699 Solingen, Germany
Phone : +49 (0) 212 232380 Fax : +49 (0) 212 2323899

DermaI Curette Biopsy Technique for Various Dermatological Diseases

Manabu MAEDA, MD & Ph D

Department of Dermatology, Prefectural Gifu Hospital, 4-6-1 Noishiki, Gifu City, 500-8717, Japan

In diagnosing various dermatological diseases, skin biopsies are not just essential for making a diagnosis and determining whether there is a malignant tumor, but they are the techniques^{1,2)} that are required in daily clinical practice by dermatologists and doctors of other departments. However, the common and frequently performed biopsies using a surgical knife sometimes require imperative experience in determining incision directions and suturing the skin, and therefore may not be suitable in an emergency or for unskilled non-dermatologists in various situations. In such cases, skin biopsies or minor surgery⁴⁾ by the shave removal method³⁾ using a blade or those using a dermal curette like a biopsy punch would be relatively easy and useful.

This paper reports in detail biopsy techniques using a disposable dermal curette and introduces clinical applications of those techniques.

[Methods and Techniques]

· List of required devices

Adson type hooked forceps, camper scissors for ophthalmologic use, a syringe for local anesthesia, 5 kinds of disposable dermal curettes (Dermal Curette, Kai Industries Co., Ltd.) with diameters of 2, 3, 4, 5 and 7 mm, a diamond needle holder, and a suture thread (preferably on a needle).

1) Skin biopsy

(to determine whether there is a malignant tumor)

First, we wanted to determine of the desired site for the skin biopsy (arrow in Fig. 1a). Local anesthesia is applied to this cross-marked site using lidocaine containing 1% epinephrine. It is, however, important to intracutaneously inject the anesthetic lidocaine into a shallow depth in the peripheral part or bank around the center rather than the central part (Fig. 1b). Usually an amount of about 1 to 2 cc is appropriate, but the amount of local anesthetic lidocaine varies depending on the patient, while an additional injection will be necessary if the appropriate effect is not achieved. In order to prevent necrosis, it is beneficial to avoid using epinephrine at distal sites such as fingertips and toe tips. After the effect of anesthesia has been ascertained with hooked forceps, a skin specimen is collected by shaving the skin surface with the disposable dermal curette (Fig. 1c). At this time, the excessive pressing of the curette might crush or deform the biopsy specimen. The key to a successful biopsy is to pull the blade toward you like you're peeling a potato or a cucumber, while giving proper force. The skin specimen is quickly collected from the depth of dermis using the above method. It is necessary to pay close attention to see whether there are any important organs, blood vessels, or nerves running into the inferior bed. After the skin specimen has

been collected, the wound may be sutured as shown in Fig. 1d. In some cases, the progress of the open wound may be observed as it is, by merely applying antibiotics containing ointment. The use of a suture thread (e.g., Ethylon®) on a needle would be more convenient as shown in the figure. In general, it is ideal that the muscular orientation be vertical to the incision and suture line, however this principle is not always applicable depending on the site. Therefore, the direction in which the skin can be easily gathered should be chosen to suture the wound. This method is recommended, especially in a site where tension is applied.

2) Removing skin necrosis

When various kinds of necrotic lesions such as diabetic ulcers and low-temperature burns are excised, i.e., debridman is performed; they are typically removed little by little using hooked forceps, a cooper, or a surgical knife. However, the use of this dermal curette would allow this removal to be carried out with ease and accuracy, thus shortening the time required. Of course, the site to be excised must be locally anesthetized in advance with 1% lidocaine. Fig. 2 shows a case of diabetic ulcers, where conduction anesthesia was performed and then dirty granulation was sufficiently scraped with a curette. For further details, please refer to the figure. In this case, artificial dermis (Pelnak®) was used immediately after debridman, a thin-layer skin graft was performed 2 weeks later, and then the ulcers were roughly cured 4 weeks thereafter.


Fig.1. Skin biopsy for diagnosing nevus pigmentosus and determining the degree of its malignancy: the area around the black protruded small node (indicated with an arrow in a) located almost at the central part was locally anesthetized with 1% lidocaine (b), the neck of the small node was quickly separated for a biopsy with a curette of 7 mm in diameter (c), and then the wound was sutured with a suture thread on a needle (d).


Fig.2. Diabetic ulcers developed at the right 1st toe tip: Since ulcers were covered with dirty granulation, they were difficult to epithelialize. The surrounding epidermis that thickened like a bank is shaved with a curette of 2 mm in diameter after conduction anesthesia was performed. Subsequently, an attempt was made to excise the poorly proliferated granulation at the central part with a curette of 3 mm in diameter (c). The size of the ulcerous area was traced with a piece of sterilized filter paper. Artificial dermis that was cut out to that precise size was adhered to the filter paper, sutured, and then covered by tie-over dressing (d). The tie-over dressing was removed 5 days later, an ultra thin-layer skin graft was performed 2 weeks later, and then the wound was almost cured 4 weeks following (the figure is omitted).


Fig.3. Treatment for plantar clavus: The blade of a curette with a diameter of 2 mm is positioned at the part equivalent to the eye of the clavus (b), and then it's rotated 1 full turn to shave the central part (c). Immediately after this, salicylic acid (Speel®) plaster is applied to only the central part. Several days later, the whitened macerated site is similarly treated with the curette.


Fig.4. Treatment for clavus with a cutaneous horn: A clavus, extruding like a cutaneous horn, can be treated similarly as shown in Fig. 3. The cutaneous horn part is pressed against the blade of a curette and then the curette is rotated 1 full turn to shave the central part. This allows for the obtaining of satisfactory results as follows; A, before treatment for the plantar clavus; B and C; the curette is being used; a; clavus in toe before treatment, and lastly b; after treatment.

3) Treatment for clavus and tylosis

To treat a tylosis, a sufficiently satisfactory result can be obtained by merely shaving it little by little, like peeling a cucumber with a curette. For a plantar clavus, however, as shown in Fig. 3a, more skill is required. Shaving the central eye of the clavus with a curette (Figs. 3b, c, d) opens a round aperture. Salicylic acid plaster is again applied to this aperture by letting the plaster flow into itself. The lesion is excised again by the same technique several days later. This method is safe and effective. In particular, this technique is very effective for a clavus that has high keratinization and protrudes such as a cutaneous horn (Fig. 4A, a). For a clavus that involves a secondary infection, a satisfactory result may not be obtained and this method should be avoided here whenever possible.

Discussion and Conclusions

As stated above, skin biopsies are one of the more essential techniques for dermatological diseases. In daily dermatological practice, skin biopsies must be performed accurately in a short period of time, with method that meets this requirement as naturally being regarded as the best. Skin biopsies using an ordinary surgical knife have a problem of requiring time and experience ^{1,2)}. However, this method has an advantage of only requiring a syringe for local anesthesia, a dermal curette, hooked forceps, a diamond needle holder for suture, and a suture thread without requiring any conventional surgical knife. Presently, the above method is frequently used as a simple method in daily dermatological practice. It is applicable especially for superficial lesions. Dermal curettes with a variety of diameters

including 2, 3, 4, 5, and 7 mm have become commercially available, making this medical device applicable in a very broad range of situations when the specimen collection site and the patient's age are taken into consideration. It must be emphasized that this method can be easily applied to not only the incision of a clavus or a necrotic lesion that requires treatment within a very limited time in daily dermatological practice, but to that of a tylosis. Skin biopsies and minor surgeries using a dermal curette together with their brief descriptions have been introduced above. In time, they will be of great use in daily dermatological practice.

References

- 1) Alguire PC, Mathes BM. Skin biopsy techniques for the internist. J General Int Med. 13; 46-54, 1998
- 2) Todd P, Garioch JJ, Humphreys S, Seywright M, Thompson J, du Vivier AW. Evaluation of the 2-mm punch biopsy In dermatologic diagnosis. Clin Exp Dermatol. 21; 11-13, 1996
- 3) Harvey DT, Fenske NA. The razor blade biopsy technique. Introduction of the adaptor-designed shave biopsy instrument. Dermatol Surg. 21;345-7, 1995
- 4) Zalla MJ. Basic cutaneous surgery. Cutis. 54; 172-186,1994

DISPOSABLE DERMAL CURETTE
sterilized

- Provided with a stainless steel blade that is suitable for shaving soft tissues in the affected area.
- Disposable ensuring that reuse involving a risk of infection can be avoided.

- Dermatological treatment
- Dermatological tissue test
- Animal experiments

CE
0197

Medical device permission of Japan No. 21BZ0070

Manufactured by

kai industries co.,Ltd.
1110, Oyana, Seki City, Gifu Prefecture, 501-3992 Japan
Phone : +81 (0) 575 28 6600
Fax : +81 (0) 575 28 6611
Email : info-ind@kai-group.com
http : //www.kai-ind.co.jp

Punch (Trepan) Biopsy Technique for Various Dermatological Diseases

Manabu MAEDA, MD & Ph D

Department of Dermatology, Prefectural Gifu Hospital, 4-6-1 Noishiki, Gifu City, 500-8717, Japan

In diagnosing various dermatological diseases, skin biopsy is an essential technique for making a diagnosis in daily clinical practice by dermatologists and doctors of other departments. Since this common biopsy is performed with a surgical knife, and only after essential experience is gained in determining incision directions and suturing of the skin, it may not be suitable for an emergency situation or for unskilled non-dermatologists for various reasons. In these such cases, the shave removal method biopsy using a blade or biopsy punch would be relatively easy and useful. This paper reports on the detailed biopsy techniques using a disposable biopsy punch and introduces the kind of tips of biopsy punches with application methods in daily dermatological practice.

[Methods and Techniques]

• List of required devices

Adson type hooked forceps, camper scissors for ophthalmologic use, a syringe for local anesthesia, a total of 10 kinds of disposable biopsy punches (Biopsy Punch, KAI Industries Co, Ltd.) with diameters of 1, 1.5, 2, 2.5, 3, 3.5, 4, 5, 6, and 8 mm, a diamond needle holder, and a suture thread (preferably on a needle).

1) Skin biopsy

(to determine whether there is a malignancy or not)

A Delmark® (marker exclusively for dermatological use) is used to mark a cross at the center of the desired site for a skin biopsy (Fig. 1b). Local anesthesia is applied to this cross-marked site using 1% lidocaine (Xylocaine®) containing epinephrine (adrenaline). It is, however, important to intracutaneously inject the anesthetic into a shallow depth in the peripheral part or bank around the center rather than the central part (Fig. 2a). Usually an amount of about 2 cc is appropriate, but the amount of local anesthetic varies depending on the patient, while finally an additional injection will be necessary if the appropriate effect is not achieved. In order to prevent necrosis, it is beneficial to avoid using epinephrine at distal sites such as fingertips and toe tips. After the effect of anesthesia has been ascertained with the hooked forceps, the leading end of the punch is positioned at the center of the cross mark and then the punch is quickly pressed vertically into the skin surface until it reaches the dermis. At this time, by excessive pressing the punch might crush or deform the biopsy specimen. The key to a successful biopsy is to uniformly distribute force between pressing and drilling. This approach can be more easily understood if you try to imagine how to drill a hole in a board with a gimlet. It is essential that the punch be inserted until it reaches the dermis, in particular, the subcutaneous fat tissues in a spiral state. Even if the punch is inserted until the border between its plastic part and its metal blade, it is not necessary to be concerned about damaging subcutaneous tissues as long as no important organs, blood vessels, or nerves run into the inferior bed, although this varies depending on the site. If the punch completely cuts subcutaneous fat tissues at one time, they are often discharged naturally. This facilitates obtaining the cut specimen as it is raised. If this should be impossible, the leading end of the specimen must be clamped very lightly with the hooked forceps, slightly drawn out upward, and then severed from the inferior bed using the ophthalmologic camper scissors. After a successful hemostasis has been ascertained, the wound is sutured. As shown in Fig. 2, the use of a suture thread (e.g., Ethylon®) on a needle would be more convenient. In general, it is ideal if the muscular orientation is vertical to the incision and suture line, but this principle is not always applicable depending on the site. Therefore, the wound should be sutured circularly. If the wound is deformed elliptically, it should be sutured assuming that this part has been excised with a surgical knife. This method is recommended especially in a site where tension is applied.

Fig. 1. A clinical image (a) of verruca senilis (senile wart) developed on the face and a view (b) positioning of the punch biopsy site.

2) Total enucleation of small skin tumor

A relatively small lentigo or nevus cell nevus that is covered by the inner diameter of the punch can be completely enucleated with accuracy and efficiency by using a punch than by using an ordinary surgical knife. A Delmark® is used to mark a cross at the center of the site to be excised and then a punch is inserted into exactly the center. This enables perfect total enucleation without a mistake. The technique is exactly the same as described in the above 1); the details and figures are omitted.

3) Simplified removal of a subcutaneous tumor (boring method)

The central part of an epidermal cyst is bored with the punch and then by both the contents of the epidermal cyst, and subsequently the constituent walls are pushed out with fingers^{4,5}. This method is applicable for epidermal cysts that are movable without inflammation and adhesion to the inferior bed.

For epidermal cysts that involve a secondary infection or are recurrent, a satisfactory result may not be obtained and the boring method should be avoided wherever possible. However, if pus discharge is intended by incising an infectious epidermal cyst, a punch may be used because a round wound can be more easily filled with a tampon rather than a wound that has been linearly incised with a knife. This technique has a problem of requiring some essential experience and skill. However, its advantages of minimizing bleeding and the surgical wound are noteworthy. If this method is further developed, subcutaneous nodes that have no wall unlike epidermal cysts and benign tumors that exist subcutaneously can be almost completely excised. Fig. 3 shows a case where a calcifying epithelioma was almost completely removed. An even larger one could be excised using a punch of 5 to 6 mm in diameter. Since the scar was small, the cosmetic result was evaluated significantly well.

Discussion and Conclusions

As stated above, skin biopsies are one of the more essential techniques for dermatological diseases. In daily dermatological practice, skin biopsies must be performed accurately in a short period of time with the method that meets this requirement as naturally being regarded as the best. Skin biopsies using an ordinary surgical knife have a problem of requiring a person with experience and skill^{1,2}). However, this method has an advantage of only requiring a syringe for local anesthesia, a punch, hooked forceps, a diamond needle holder for suture, and a suture thread without requiring any conventional surgical knife. Presently, the above method is frequently used as a simple method in daily dermatological practice. It cannot be affirmed, however, that this method is applicable in all cases. Nevertheless, punch with a variety of diameters including 1, 1.5, 2, 2.5, 3, 3.5, 4, 5, 6, and 8 mm have become commercially available making this medical device applicable in a very broad range situations when the specimen collection site and the patient's age are taken into consideration. Thus, this method enables not just skin biopsies for infants, but even total enucleation of a plantar lentigo or nevus cell nevus that falls within the inner diameter of the punch used, as mentioned earlier. It must be emphasized that this method was reported as a simplified one for excising an epidermal cyst³) and it is applicable for the frequently used shave biopsy. A prototype of a long type punch with a longer effective neck was also used⁶). Using this allows accurately collecting a specimen during biopsy for a nodule or a lesion that is extended to the deep subcutaneous area, e.g., lupus


Fig. 2. Skin biopsy using a punch for verruca senilis (same case as in Fig. 1): (a) local anesthetic is applied to the surrounding area, (b) the punch is vertically inserted from the positioned site, (c) the punch is removed, (d) the biopsy specimen is drawn out from the inferior bed using hooked forceps and cut off by camper scissors for ophthalmologic use, (e) the suture thread with a needle is taken out, (f) the aperture created by the punch is sutured, (g) suture should be performed quite extensively from the outside, and (h) gauze is applied. Now the biopsy has been completed. Since it is critical to pathophysiologically determine whether there is a malignant tumor, a skin biopsy was performed on the central part using a punch.

erythematosus profundus or morphea. The use of a punch enables to minimize the extent of a biopsy and the scar after the biopsy, thus being regarded as a useful tool for good cosmetic results. Skin biopsies and minor surgical techniques using a punch together with their brief descriptions have been introduced above. In daily dermatological practice, more diversified methods would be made possible as they are devised. Therefore, skin biopsies using a punch will increasingly become useful.


Fig. 3. Technique for excising a calcifying epithelioma by applying the punch method matching like (f), an aperture is opened with the punch at the central part of the tumor. The leading end of the tumor is drawn out with hooked forceps or a skin hook and then excised from the surrounding tissues using ophthalmologic camper scissors. This is a useful technique to know and retain.

References

- 1) Alguire PC, Mathes BM. Skin biopsy techniques for the internist. *J General Int Med.* 13; 46-54, 1998
- 2) Todd P, Garioch JJ, Humphreys S, Seywright M, Thompson J, du Vivier AW. Evaluation of the 2-mm punch biopsy in dermatologic diagnosis. *Clin Exp Dermatol.* 21; 11-13, 1996
- 3) Harvey DT, Fenske NA. The razor blade biopsy technique. *Introduction of the adaptor-designed shave biopsy instrument.* *Dermatol Surg.* 21;345-7, 1995
- 4) Kamide R, Atheroma; *Dermatological Practice 4, Practice for Day Surgery* (Edited by Ohara K, Miyaji Y, Takigawa M), Bunkodo, Tokyo, 1998, p183-187
- 5) Maeda M, Fujisawa T, Araki M, Yamazaki T, Simplified operation method using a punch instrument for epidermal cysts. XXIII Congress of the international society for dermatologic surgery and national meeting of Portuguese group of dermatologic surgery. Sept.18-22, 2002, at Porto
- 6) Maeda M, Satoh M, Iwata H, Yamazaki T, Sawada Y, Araki M, Clinical application and usage of long type disposable biopsy punch (KAI). *Nishinihon J Dermatol* 62; 783-787, 2000

DISPOSABLE BIOPSY PUNCH

sterilized

- Available in a broad range of sizes.
- Disposable ensuring that reuse involving a risk of infection can be avoided.

- Dermatological treatment
- Dermatological tissue test
- Animal experiments

● The long type is available in a blade length of 15 mm (standard type: 7 mm).

CE
0197

Medical device permission of Japan No. 21BZ0070

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Riferimenti

- 1) Algurs PC, Mathes BM Skin biopsy techniques for the internist. J General Int Med. 13; 46-54, 1998
- 2) Todd P, Garioch JJ, Humphreys S, Seywright M, Thompson J, du Vivier AW. Evaluation of the 2-mm punch biopsy in dermatologic diagnosis. Clin Exp Dermatol. 21; 11-13, 1996
- 3) Harvey DT, Fenske NA. The razor blade biopsy technique. Introduction of the adaptor-designed shave biopsy instrument. Dermatol Surg. 21; 345-7, 1995
- 4) Kamide R, Atheroma; Derrnamological Practice 4, Practice for Day Surgery (Edited by Ohara K, Miyaji Y, Takigawa M), Bunkodo, Tokyo, 1998, p. 183-187
- 5) Maeda M, Fujisawa T, Araki M, Yamazaki T. Simplified operation method using a punch instrument for epidermal cysts. XXIII Congress of the international society for dermatologic surgery and national meeting of Portuguese group of dermatologic surgery. Sept. 18-22, 2002, at Porto
- 6) Maeda M, Satoh M, Iwata H, Yamazaki T, Sawada Y, Araki M, Clinical application and usage of long type disposable biopsy punch (KAI), Nishinohon J Dermatol 62; 783-787, 2000

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